

# H6 plus

## Closed Urine Drainage System



### Product features

In order to keep the risk of nosocomial infections to a minimum HSO offers a **closed sterile urine drainage system** – the H6 plus.

The H6 plus urine drainage system is easy to use and in addition to the normal standards also offers a number of user-friendly extras. An additional blue clamp on the transition tube helps fixation. The bottom outlet valve opens and closes exactly and is drip-free. Needleless taking of urine samples. The complete connection piece to the catheter, including emptying port, is easy to clean hygienically.

The ventilated urine bag has a capacity of 2,000 ml and can be fixed in place with either the cord or hook. The bacteria-proof diaphragm drip chamber prevents retrograde bacteria migration and the risk of infection that this involves. The overall design of the drip chamber also prevents urine from running back into the catheter.

The entire tube is kink-proof. The system is available in two versions: with a hard and with a soft tube.

## Instructions for use



The H6 plus urine drainage system has a **hygienic port for needleless taking of samples** and a **drip-free bottom outlet valve**.

## Benefits of the H6 plus at a glance

- Needleless taking of urine samples
- Bacteria-proof drip chamber with diaphragm
- Exact and drip-free bottom outlet valve
- Kink-proof transition tube (120 cm) with slide and fixation clamp
- Fixation with either hook or cord
- Ventilated urine bag (2,000 ml)
- Easy-to-read markings
- Latex-free



## Accessories

### H3 Urine Collecting Bag

For safe and hygienic emptying of urinary drainage systems as well as for urine measurement and urine disposal.

- Very flexible funnel connector – optimally adjusts to the anti-drip bottom valve
- Sliding clamp on the tube for opening and closing the urine collecting bag
- Four eyelets for fastening/hanging up
- 2,500 ml capacity
- Disposable product, latex-free, non-sterile

### ORDER OPTIONS HSO URINE DRAINAGE SYSTEM H6 PLUS

Product	Shipping Unit	Order No.	
Closed urine drainage system with hard tube	25 pieces (single packed, sterile)	HEICSB2000P	CE 0408
Closed urine drainage system with soft tube	25 pieces (single packed, sterile)	HEICSB2000P-S	CE 0408
H3 Urine Collecting Bag	100 pieces	HSOH3AB2500U	CE