

Over-bed table

Over-bed table with inclinable table tray



Product features

With its **continuous height adjustment** and **functional table trays** the **HSO over-bed table** always offers the best solution for both patients and personnel alike. In addition to its **practical design** the HSO over-bed table also has an inclinable table tray with an **anti-slip bar** for even more security.

The static table tray has two handles that allow **simple manoeuvring** of the table. The rubber lined casters – two of which are lockable – ensure sufficient stability.

The HSO over-bed table is the **optimum solution** for use in the intensive care unit, **especially in the recovery room**, and the **ideal accessory in the patient's room**.

Technical Data

Table tray	
Inclinable table tray	640 x 400 mm
Inclination angle	3 settings: 30°/45°/60° (+/- 2°)
Static table tray	200 x 400 mm
Table tray thickness	8 mm



Design	
Height adjustment	Continuous, 790-1,035 mm, with damping mechanism
Table column design	Column and chassis connected with a welded joint & nut/bolt
Finish	Powder coated finish
Casters	4 casters each with a 50 mm diameter, 2 of which are lockable



ORDER OPTIONS HSO OVER-BED TABLE		
Product	Shipping unit	Order No.
Over-bed table	1 piece	HS0355-685-036

